

Sealver GEL

SENSITIVE PROFESSIONAL REPAIR



NO MORE SENSITIVE TEETH



WATCH THE VIDEO
FOR PROPER
USE

 **DOCTOR
SMILE**
doctor-smile.com

SAY GOODBYE TO DENTINAL HYPERSENSITIVITY



Sealver GEL can make your patients forget this disease in just **3 minutes**: thanks to its localized effect, the gel activated by the laser light desensitizes and seals the dentinal tubes liable for this pathology while protecting the glaze and the dentin. **The benefit coming from Sealver GEL is immediate**, it spreads positive feels to the patient thanks to a quick treatment: a premium service which can make the difference also in the subsequent stages of the dental bleaching.

FORGET THE DISCOMFORT
OF HEAT AND COLD



FREE FROM DENTINAL HYPERSENSITIVITY
IN JUST 3 MINUTES



QUICK



NO PROBLEM



WOW EFFECT



PROFESSIONAL



LASTING

 **DOCTOR
SMILE**



LAMBDA SpA

Via dell'Impresa | 36040 Brendola (VI) | Italy
T +39 0444 349165 | F +39 0444 349954
info@lambdaspa.com | lambdaspa.com



MANI MEDICAL GERMANY GmbH
Hertha-Sponer-Straße 2
61191 Rosbach v. d. Höhe | Germany



MODEL DEVICE

INFO

· doctor-smile.com
· info@doctor-smile.com

